

DISPOSABLE FLOVAC® LINER, 1000 ML FOR ASPIRET, ASKIR 20, 30



Ref. BCASPI99



25/08/2024 20:38:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DISPOSABLE FLOVAC® LINER, 1000 ML FOR ASPIRET, ASKIR 20, 30



Ref. BCASPI99

Description

- Flexible disposable polyethylene pouch
 - Hermetically welded to a rigid lid inside which a hydrophobic and antibacterial anti-reflux filter is placed. This filter also works as an overflow, suction valve deactivation system whenever the maximum filling capacity is reached.
 - The pocket must be placed in its FLOVAC® rigid polycarbonate autoclavable CONTAINER.
- FLOVAC® systems therefore replace the use of autoclavable jars, which reduces cleaning and sterilization costs.
- Compatible with New Aspiret and Series Askir Suction Units. (Not compatible with C30 series)

Features

CE Marking	CE0123	Medical device	Class Im (measuring function)
Bar code	3661040151137	HS Code	90189084
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.2
Bare Unit / Weight (kg)	0.2		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:38:03 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr