

## SPO2 SENSOR PAEDIATRIC CLIP FOR COMEN NC3/NC5/C31

# COMEN

Ref. DBCPEDI01



CE  
0537

DISPOSITIF MEDICAL  
CLASSE  
**IIb**  
MEDICAL DEVICE

GARANTIE  
1 YEAR  
**1**  
1 YEAR  
GUARANTEE

 25/08/2024 21:07:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

# SPO2 SENSOR PAEDIATRIC CLIP FOR COMEN NC3/NC5/C31

# COMEN

Ref. DBCPEDI01

Designed for accurate and fast measurement of oxygen saturation

## Description

Technical specifications: % oxygen 90-100% accuracy +/- 2% 70-90%

accuracy +/-3%  
<70% not specified

## Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219424	HS Code	90181990
Warranty (months)	12	Country	Uruguay
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 25/08/2024 21:07:17 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)