

SPO2 SILICON WRAP SENSOR FOR COMEN NC3/NC5/C34

COMEN

Ref. DBSWRAP01



25/08/2024 21:08:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SPO2 SILICON WRAP SENSOR FOR COMEN NC3/NC5/C34

COMEN

Ref. DBSWRAP01

Designed to ensure accurate and fast measurement of oxygen saturation.

Description

Technical specifications % oxygen 90-100% Accuracy +/-2% 70-

90% Accuracy +/-3%
<70% not specified

Features

CE Marking	CE0537	Medical device	Class IIb
Bar code	3661040219455	HS Code	90181990
Warranty (months)	12	Country	Uruguay
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	13	Packed unit / Height (cm)	3
Packed unit / Length (cm)	13	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 25/08/2024 21:08:29 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr