

3 INCHES DISPOSABLE PADS FOR TDH39 AND DD45 HEADPHONES (X500)


Ref. DSIPROT01



 25/08/2024 21:09:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

3 INCHES DISPOSABLE PADS FOR TDH39 AND DD45 HEADPHONES (X500)

Ref. DSIPROT01

Description

Disposable protectors for helmets. They are compatible with DD45 and THD-39 helmets and protect the helmet between each patient for optimal hygiene.

Packaging: box of 500

Features

CE Marking	NONE	Medical device	No
HS Code	42032910	Warranty (months)	12
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	30	Packed unit / Height (cm)	15
Packed unit / Length (cm)	30	Packed unit / Weight	0.8
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 25/08/2024 21:09:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr