


RECHARGEABLE BATTERY FOR MCC/MPE/MPD/MPC/MBC

KERN[®]


Ref. EKMB08



 25/08/2024 20:42:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

RECHARGEABLE BATTERY FOR MCC/MPE/MPD/MPC/MBC



Ref. **EKMBC08**

Description

Operation with internal battery, service life up to 40 h without backlight, charging time approx. 12 h Suitable for KERN MWA wheelchair weighing platform, KERN MTA scales with retaining bar, KERN MPE, KERN MPD and KERN MPC column scales, KERN MCC scale and KERN MBC baby scales

EKMBC08 suitable for serial numbers WF1434431 and higher (otherwise see Battery EKMBC05)

Features

CE Marking	NONE	Medical device	No
Bar code	3661040173535	HS Code	85075000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.184
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:42:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr