


ADULT TIPS FOR PEAK FLOWMETERS Ø 30 MM - BOX OF 100 PCS



Ref. L35587

CE



 25/08/2024 21:02:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ADULT TIPS FOR PEAK FLOWMETERS Ø 30 MM - BOX OF 100 PCS



Ref. **L35587**

Description

Adult tips for SPENGLER peak flow meter.
 Dimensions of the tips 30mm.
 Sold in boxes of 100 pieces.

Features

CE Marking	CE	Medical device	Class I
Bar code	3700446001681	HS Code	48229000
Country	China		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.83
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 21:02:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr