

S/S BOX 28 X 14 X 6 CM

Ref. AI28146



25/08/2024 19:09:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S/S BOX 28 X 14 X 6 CM



Ref. AI28146

Convenient for storing surgical instruments

Description

- Sold with lid to protect instruments
- Decontamination with detergent and disinfectant -
- Sterilizable by autoclave (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

Features

| | | | |
|------------------------------------|---------------|-------------------------|----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040019574 | HS Code | 90189090 |
| Warranty (months) | 12 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 14 | Bare Unit / Width (cm) | 14 |
| Packed unit / Height (cm) | 6 | Bare Unit / Height (cm) | 6 |
| Packed unit / Length (cm) | 28 | Bare Unit / Length (cm) | 28 |
| Packed unit / Weight | 0.872 | Bare Unit / Weight (kg) | 0.872 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

25/08/2024 19:09:14 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr