

KOCHER FORCEPS 1X2 TEETH STRAIGHT 13 CM

Ref. IP02013



25/08/2024 20:00:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

KOCHER FORCEPS 1X2 TEETH STRAIGHT 13 CM



Ref. IP02013

Allows you to clamber or fix elements (tubing, etc.). Used mainly on sturdy fabrics.

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021379	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Bare Unit / Width (cm)	9
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.3
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	19.5
Packed unit / Weight	0.03	Bare Unit / Weight (kg)	0.029
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:00:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr