

## KOCHER FORCEPS 1X2 TEETH STRAIGHT 16 CM

Ref. IP02016



25/08/2024 20:01:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## KOCHER FORCEPS 1X2 TEETH STRAIGHT 16 CM



Ref. IP02016

Allows you to clasper or fix elements (tubing, etc.). Used mainly on sturdy fabrics.

### Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021409	HS Code	90189090
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	18	Bare Unit / Width (cm)	7.2
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.4
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	16.5
Packed unit / Weight	0.05	Bare Unit / Weight (kg)	0.044
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

25/08/2024 20:01:10 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)