

ROCHESTER-PEAN PLIERS - KOCHER, S/G, RIGHT, 24 CM

Ref. IP02224



25/08/2024 20:04:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ROCHESTER-PEAN PLIERS - KOCHER, S/G, RIGHT, 24 CM



Ref. IP02224

Convenient to stop bleeding without damaging fragile tissues during an operation. Also used to perform care (dressings, etc.).

Description

Holtex instruments are all made of top-quality stainless steel (AISI 304, 410, 420 and 420A, depending on the type of instrument) to guarantee good durability and resistance to corrosion.

The steels come into contact with human tissues and must therefore comply with international standards and regulations in force.

Holtex guarantees that its medical devices comply with the 2017/745 EU regulation and ISO 13485.

100% IN-STORE reference: this reference of service instrumentation is guaranteed to be available immediately.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040021614	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	28	Bare Unit / Width (cm)	8.1
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.5
Packed unit / Length (cm)	9	Bare Unit / Length (cm)	24.2
Packed unit / Weight	0.08	Bare Unit / Weight (kg)	0.083
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

25/08/2024 20:04:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr