

KUNSCHER NAIL EXTRACTOR (WITHOUT HOOKS)



Ref. IC01416



26/08/2024 20:50:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

KUNSCHER NAIL EXTRACTOR (WITHOUT HOOKS)



Ref. IC01416

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040008707	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	6	Packed unit / Height (cm)	2
Packed unit / Length (cm)	20	Packed unit / Weight	0.35
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 20:50:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr