

SURGERY SET WITH ALUMINIUM BOX 18X9X3 CM



Ref. IA05601



26/08/2024 19:31:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SURGERY SET WITH ALUMINIUM BOX 18X9X3 CM



Ref. IA05601

Description

Consisting of:

- IC09411 Straight iridectomy scissors, 11.5 cm
- IC11414 Straight round-ended scissors, 14 cm
- ID05403 no.3 scalpel handle
- IP02014 Straight kocher toothed forceps, 14 cm
- IP21414 Pean Murphy forceps, 14 cm
- IP26014 toothed dissection forceps, 14 cm
- IP26014 Dissection forceps without teeth, 14 cm
- IP26410 Splinter forceps, 10 cm
- IP29414 Mayo-Hegar needle holder, 14 cm
- IS09414 Cotton holding stylet, 14 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040006314	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	9
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	18
Packed unit / Weight	0.45	Bare Unit / Weight (kg)	0.45
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 19:31:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr