

REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 4 MM), 15CM



Ref. IA08412



26/08/2024 19:35:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

REDON INTRODUCING NEEDLE STRAIGHT (DIAM. 4 MM), 15CM



Ref. IA08412

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040006420	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	11.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	0.3
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	16.6
Packed unit / Weight	0.022	Bare Unit / Weight (kg)	0.022
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 19:35:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr