


ALUMINIUM BOX 18X9X4 CM, GREY

Ref. AAT1894



26/08/2024 17:58:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ALUMINIUM BOX 18X9X4 CM, GREY



Ref. AAT1894

Convenient for storing surgical instruments

Description

- Hot sterilization or autoclave (134°-18 minutes)
- Do not wash in an ultrasonic cleaner or in a washing machine
- Warning: aluminum is incompatible with some decontamination products

Format 18 x 9 x 4 cm

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040000275	HS Code	90184990
Warranty (months)	12	Country	France
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	9	Bare Unit / Width (cm)	9.3
Packed unit / Height (cm)	2	Bare Unit / Height (cm)	2.5
Packed unit / Length (cm)	18	Bare Unit / Length (cm)	18.3
Packed unit / Weight	0.157	Bare Unit / Weight (kg)	0.157
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 17:58:26 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr