

S/S KIDNEY DISH 16 X 8.5 X 3.5CM

Ref. AI16124



26/08/2024 18:28:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S/S KIDNEY DISH 16 X 8.5 X 3.5CM



Ref. AI16124

Ideal for storing instruments or temporarily installing medical waste

Description

- Decontamination with detergent and disinfectant - Autoclavable (134°-18 minutes)- In the " ATTACHED DOCUMENTS" tab, you can download the "Cleaning and sterilization" manual
- Reference 100% IN-STORE: this service instrumentation reference is guaranteed available immediately
- Available in different sizes

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019468	HS Code	90189090
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Bare Unit / Width (cm)	8.5
Packed unit / Height (cm)	7	Bare Unit / Height (cm)	3.5
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	16
Packed unit / Weight	0.077	Bare Unit / Weight (kg)	0.077
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 18:28:07 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr