

S/S BOX PERFORATED 20 X 10 X 5CM PERFORATED TOP / BOTTOM

Ref. AI20105P



26/08/2024 18:32:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

S/S BOX PERFORATED 20 X 10 X 5CM PERFORATED TOP / BOTTOM



Ref. AI20105P

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067742	HS Code	90189090
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	12	Bare Unit / Width (cm)	10.5
Packed unit / Height (cm)	6	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	20.7
Packed unit / Weight	0.407	Bare Unit / Weight (kg)	0.407
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

📅 26/08/2024 18:32:19 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

🌐 www.gsh-med.fr

☎ +33(0)4 42 90 31 31

✉ info@gsh-med.fr