

## S/S BOX PERFORATED 35 X 18 X 8 CM PERFORATED TOP / BOTTOM

Ref. AI35188P



26/08/2024 18:37:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## S/S BOX PERFORATED 35 X 18 X 8 CM PERFORATED TOP / BOTTOM



Ref. AI35188P

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040067872	HS Code	90189090
Warranty (months)	12	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	19	Bare Unit / Width (cm)	18
Packed unit / Height (cm)	9	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	38	Bare Unit / Length (cm)	36
Packed unit / Weight	1.142	Bare Unit / Weight (kg)	1.142
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

26/08/2024 18:37:56 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)