

LUBRICATED DISPOSABLE SHEATHS FOR ELEC. THERMOMETERS, PACK OF 1000

Ref. [CT40251](#)

CE



 26/08/2024 18:58:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LUBRICATED DISPOSABLE SHEATHS FOR ELEC. THERMOMETERS, PACK OF 1000

Ref. **CT40251**

Description

- Disposable case, double protection, for hygiene when taking temperature.
- The inner case prevents the thermometer from coming into contact with the mucosa.
- The color part facilitates the opening of the inner case to place the thermometer.
- The thin and resistant film guarantees accurate measurement results.
- The lubricant on the inner case facilitates the rectal introduction of the thermometer.
- For single use, the outer case protects nursing staff and worktops from contamination.
- Duration of the case: 3 years.

Available in non-lubricated version.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040069081	HS Code	9025900090
Country	Sweden		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	15	Bare Unit / Width (cm)	15.5
Packed unit / Height (cm)	8	Bare Unit / Height (cm)	8
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	21
Packed unit / Weight	0.841	Bare Unit / Weight (kg)	0.841
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 26/08/2024 18:58:09 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr