

## S/S BED PAN WITHOUT COVER SLIPPER TYPE 48 X 30CM



Ref. AI14830



26/08/2024 18:15:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## S/S BED PAN WITHOUT COVER SLIPPER TYPE 48 X 30CM



Ref. AI14830

### Description

Capacity 6 Liters

Maximum weight supported: 90kg

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040019352	HS Code	90189090
Warranty (months)	12	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	55	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	8	Packed unit / Length (cm)	30
Bare Unit / Length (cm)	48	Packed unit / Weight	1.08
Bare Unit / Weight (kg)	1.07		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

26/08/2024 18:15:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)