

LID FOR BEDPAN (STAINLESS STEEL) 48X30CM



Ref. AI14831



26/08/2024 18:15:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

LID FOR BEDPAN (STAINLESS STEEL) 48X30CM



Ref. AI14831

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040032917	HS Code	90189090
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	54	Bare Unit / Width (cm)	30
Packed unit / Height (cm)	2	Packed unit / Length (cm)	39
Bare Unit / Length (cm)	48	Packed unit / Weight	0.416
Bare Unit / Weight (kg)	0.416		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

26/08/2024 18:15:46 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr