

KILLIAN NASAL SPECULUM 50 MM WITH LATERAL SCREW



Ref. IS23050



27/08/2024 06:55:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

KILLIAN NASAL SPECULUM 50 MM WITH LATERAL SCREW



Ref. IS23050

Features

CE Marking	CE	Medical device	Class Ir
Bar code	3661040031156	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	5.5
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	5.2
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	15
Packed unit / Weight	0.099	Bare Unit / Weight (kg)	0.099
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 06:55:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr