

STAINLESS STEEL BED PAN SLIPPER TYPE 42X29 CM WITHOUT LID



Ref. AI14832



29/08/2024 12:28:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

STAINLESS STEEL BED PAN SLIPPER TYPE 42X29 CM WITHOUT LID



Ref. AI14832

Description

- Maximum weight supported: 90 kg

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040115207	HS Code	90189090
Warranty (months)	12	Country	Pakistan
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	30	Bare Unit / Width (cm)	29
Packed unit / Height (cm)	17	Packed unit / Length (cm)	56
Bare Unit / Length (cm)	42	Packed unit / Weight	1.08
Bare Unit / Weight (kg)	1.481		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 12:28:02 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr