

CFS-A01 INTERFACE CABLE RS-232 TO CONNECT AN EXTERNAL DEVICE




Ref. EKCFS01

 29/08/2024 11:54:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

CFS-A01 INTERFACE CABLE RS-232 TO CONNECT AN EXTERNAL DEVICE



Ref. **EKCFS01**

Description

Cable length approx. 1 m

Features

CE Marking	NONE	Medical device	No
Bar code	3661040108421	HS Code	84231090
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 11:54:48 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr