

WHEELCHAIR PLATFORM SCALE MWS

KERN®


Ref. EKMWS03



CE
0297

DISPOSITIF MEDICAL
CLASSE
Im
MEDICAL DEVICE

GARANTIE
3 ANS
3
3 YEARS
GUARANTEE



 29/08/2024 11:57:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

WHEELCHAIR PLATFORM SCALE MWS



Ref. EKMWS03

Large wheelchair weighing platform, lowered for easy access of electric wheelchairs

Description

- Particularly suitable for weighing patients in wheelchairs thanks to the low-slung platform accessible from both sides with tare weight memorization of up to five (electric) wheelchairs
- The large steel weighing pan is ideal for weighing strong, overweight patients
- Stable and non-slip thanks to its height-adjustable rubber feet
- Hold function: for agitated patients, the average value looks for a stable weighing value and "freezes"
- BMI function for calculating insufficient/normal/overweight weight
- 2 casters and 1 handle make it easy to carry and store it for a small footprint (unless barriers are added)
- LCD screen, height of digits 25 mm
- Graduation 0.1 kg
- Battery operated (6 x 1.5V AA not included) or mains (power supply supplied) or optional with rechargeable battery pack
- Homologable in Class III
- Medical approval according to 93/42/EEC

Technical specifications:

- Display dimensions : L x P x H 210 x 110 x 50 mm
- Total dimensions : L x P x H 1155 x 800 x 66 mm
- Display cable length approx. 180 cm
- Net weight approx. 28.4 kg

For any measuring device requiring medical Class III approval, the accessories must be ordered at the same time as the device so that the approval can be valid.

 29/08/2024 11:57:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

WHEELCHAIR PLATFORM SCALE MWS



Ref. **EKMWS03**

Features

CE Marking	CE0297	Medical device	Class Im (measuring function)
Bar code	4045761096707	HS Code	84231090
Warranty (months)	36	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	80
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	6.6
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	115.5
Packed unit / Weight	73	Bare Unit / Weight (kg)	42.2
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 11:57:41 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr