

HANDRAIL SET

Ref. EKMWS04

KERN[®]



 29/08/2024 11:56:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

HANDRAIL SET



Ref. **EKMWS04**

Description

- Set of adjustable retaining bars, to mount on the scale
- L 111.5cm x D 6cm x H 9 cm approximately.
- Includes: 2 retaining bars, a stabilizer stem, fixing material.
- Fits on the MWS model, with a capacity of 300kg (ref. EKMWS01)

For any measuring device requiring medical class III approval, the accessories must be ordered at the same time as it so that the approval can be valid.

Features

CE Marking	NONE	Medical device	No
Bar code	4045761096721	HS Code	84231090
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 11:56:42 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr