

## PRE-CONNECTED ELECTRODE

Ref. SJACCE01

SaverOne®



CE



27/08/2024 09:00:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## PRE-CONNECTED ELECTRODE



Ref. **SJACCE01**

### Description

For patients aged at least 8 years or weighing more than 25 kg.  
Cable length: 1.20 m

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040055145	HS Code	90181990
Country	Italy		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	17	Bare Unit / Width (cm)	18
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	1
Packed unit / Length (cm)	20	Bare Unit / Length (cm)	21
Packed unit / Weight	0.128	Bare Unit / Weight (kg)	0.128
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

27/08/2024 09:00:12 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)