

PEDIATRIC ELECTRODE

Ref. SJACCE03

Saver One®



CE



27/08/2024 09:00:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

PEDIATRIC ELECTRODE



Ref. SJACCE03

Description

Intended for children under 8 years weighing less than 25 kg.

Features

CE Marking	CE	Medical device	Class I
Bar code	3661040055169	HS Code	90181990
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	17	Bare Unit / Width (cm)	17.5
Packed unit / Height (cm)	3	Bare Unit / Height (cm)	5
Packed unit / Length (cm)	23	Bare Unit / Length (cm)	25.5
Packed unit / Weight	0.12	Bare Unit / Weight (kg)	0.12
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

27/08/2024 09:00:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr