

SEMI-AUTOMATIC DEFIBRILLATOR DSA SAVER ONE

Saver One®

Ref. SJSAVE07



CE
0051

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
6
6 YEARS
GUARANTEE

29/08/2024 11:41:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SEMI-AUTOMATIC DEFIBRILLATOR DSA SAVER ONE



Ref. **SJSAVE07**

Efficient and easy to use: reliable for all rescuers even without minimum training

Description

Effective and easy to use: reliable for all rescuers even without minimum training.

Thanks to its interactive help with its pictograms, flashing icons and clear voice message, it will guide you through all the stages of the cardio-pulmonary resuscitation for a safe rescue. Designed for public use, it authorises in 3 steps, the administration of a fast and effective treatment against sudden cardiac arrests in order to optimise 1st aid care before the arrival of the emergency services.

In semi-automatic mode, it will recommend to the operator to press the shock button only in the event of sudden cardiac arrest. The fully automatic mode is designed to deliver the therapeutic shock automatically without pressing a button or other interventions from the operator.

- Small and very light, it works with a long-life lithium battery; warranty of 5 years or 300 shocks.
- Automatic daily test and status indicator (LED light)
- Display Screen for a text message
- Graphic bar indicating the capacity of the battery
- Information button with bilingual mode
- Built-in Metronome for better support during cardiac massage
- USB connection
- Internal 1 Gb memory
- SD memory card (optional)
- Each defibrillator comes with:
 - 1 waterproof and shockproof carrying bag
 - 1 pair of pre-connected adult pads

(warranty 28 months)

 29/08/2024 11:41:52 - For information purposes only, can be subject to change

- 1st aid kit

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SEMI-AUTOMATIC DEFIBRILLATOR DSA SAVER ONE



Ref. **SJSAVE07**

Features

CE Marking	CE0051	Medical device	Class IIb
Bar code	3661040115016	HS Code	90189084
Warranty (months)	72	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	25	Bare Unit / Width (cm)	26
Packed unit / Height (cm)	16	Bare Unit / Height (cm)	15.9
Packed unit / Length (cm)	32	Bare Unit / Length (cm)	32
Packed unit / Weight	3	Bare Unit / Weight (kg)	3.27
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 11:41:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr