

DAE SAVER ONE PRO SEMI-AUTOMATIC DEFIBRILLATOR

Saver One®

Ref. SJSAVE11



CE
0051

DISPOSITIF MEDICAL
CLASSE
IIb
MEDICAL DEVICE

GARANTIE
6
6 YEARS
6 ANS
GUARANTEE



29/08/2024 11:43:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

DAE SAVER ONE PRO SEMI-AUTOMATIC DEFIBRILLATOR



Ref. SJSAVE11

Professional defibrillator 360 joules, versatile with 2 modes: manual or automated

Description

Saver One Pro semi-automatic: PROTOCOL from 50 to 360 joules - power-

- Reliable for all scenarios, it switches very easily from semi-automatic mode to manual mode if necessary

- In manual mode, it offers professionals the operational capacity of ACLS (Advanced Cardiac Life Support), allowing them to select the energy level and choose the moment for the shock, perform the synchronized cardioversion medical procedure, in ECG monitoring mode and configure access to device settings.

Moreover, if necessary, it can act as an AED providing semi-automatic operating capability, with visual and audible aid, for independent rescuers trained BLS

- Data visualization on large LCD color screen (12 x 8 cm)

- Robust, small and very light < 2 kg

- Protection rating (IP54)

- Rechargeable battery and its charging station

- Daily automatic test and status indicator (luminous led)

- Display screen for a text message

- Graph bar specifying the capacity of the battery

- Information button with bilingual mode

- Integrated metronome for better assistance to cardiac massage

- USB connection and IRDA port for printing on thermal printer

- Internal memory 1 Gb

Comes with:- 1

waterproof and shockproof carrying case- 1

pair of adult electrodes (26 months warranty)

- 1 1st aid kit (*)

- 1 complete set of wall signage - 1

SD memory card

(*) Kit Contents: compresses, 1 razor, 1 mouth-to-mouth face shield, 1 pair of sterile gloves, 1 pair of scissors

 29/08/2024 11:43:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand

13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

DAE SAVER ONE PRO SEMI-AUTOMATIC DEFIBRILLATOR



Ref. **SJSAVE11**

Features

CE Marking	CE0051	Medical device	Class IIb
Bar code	3661040128733	HS Code	90189084
Warranty (months)	72	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	4.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 11:43:52 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr