

MATRIX NEEDLE PRINTER PRINTOUT SUITABLE FOR A LONGTERM STORAGE

KERN®

Ref. EKACCE06



 29/08/2024 12:21:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

MATRIX NEEDLE PRINTER PRINTOUT SUITABLE FOR A LONGTERM STORAGE



Ref. **EKACCE06**

Description

- To print weighing values on normal paper, perfect accordingly for long-term archiving
- Needle printing
- Paper strips: width 57 mm
- Case dimensions L105×P180×H90 mm
- Please indicate in case of order your scale model so that the printer can be preconfigured accordingly at the factory
- Accessories: Paper rolls, width 57 mm (10 pieces), KERN 911-013-010

Features

CE Marking	NONE	Medical device	No
Bar code	3661040067063	HS Code	84231090
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1
Bare Unit / Weight (kg)	0.75		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 12:21:33 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr