

STATISTICS PRINTER GLP RECORD KEEPING (WITH DATE AND TIME)




Ref. EKACCE07



 29/08/2024 12:21:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

STATISTICS PRINTER GLP RECORD KEEPING (WITH DATE AND TIME)



Ref. **EKACCE07**

Description

- Silent thermoprinter with 0.75 lines/s and max. 40 characters per line.
- GLP record of adjustment report with date and time printing. Numerator. Statistical functions such as number and sum of all measured values, minimum and maximum value of the set, interval (max-min). Average value, standard deviations, available units: mg, g, kg, pcs. Languages: D, GB, F, PL.
- Case dimensions L×P×H 165×140×50 mm
- Option: Paper rolls (10 pieces), KERN YKB-A10

Features

CE Marking	NONE	Medical device	No
Bar code	3661040067070	HS Code	84231090
Country	Germany		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	1
Bare Unit / Weight (kg)	1		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 12:21:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr