


## HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°17/18



Ref. IS01417



Exemple de bougies Hégar double, inox.  
15 tailles disponibles, nous consulter.

27/08/2024 05:51:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## HÉGAR DOUBLE ENDED DILATOR STAINLESS STEEL 20 CM N°17/18



Ref. IS01417

### Description

Hégar Double Ended Dilator Stainless Steel

15 Available sizes, contact us

### Features

CE Marking	CE	Medical device	Class I
Bar code	3661040029627	HS Code	9018.9090
Warranty (months)	24	Country	Pakistan
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Bare Unit / Width (cm)	1.7
Packed unit / Height (cm)	nc	Bare Unit / Height (cm)	1.7
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	20.4
Packed unit / Weight	0.099	Bare Unit / Weight (kg)	0.099
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

27/08/2024 05:51:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)