

COLLIN VAGINAL SPECULUM 38 MM X 110 MM

Ref. IS04638



27/08/2024 06:25:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

COLLIN VAGINAL SPECULUM 38 MM X 110 MM



Ref. IS04638

Features

| | | | |
|------------------------------------|---------------|-------------------------|-----------|
| CE Marking | CE | Medical device | Class I |
| Bar code | 3661040030241 | HS Code | 9018.9090 |
| Warranty (months) | 24 | Country | Pakistan |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Bare Unit / Width (cm) | 3.8 |
| Packed unit / Height (cm) | nc | Bare Unit / Height (cm) | 6 |
| Packed unit / Length (cm) | nc | Bare Unit / Length (cm) | 11 |
| Packed unit / Weight | 0.274 | Bare Unit / Weight (kg) | 0.274 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | NA | Max temperature | NA |
| Minimum humidity | NA | Maximum humidity | NA |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

27/08/2024 06:25:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr