

## ROLL OF SILICONE TUBE Ø6X10MM LENGTH 10M FOR ASPIRET, ASKIR



Ref. BCASPI143



29/08/2024 17:35:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## ROLL OF SILICONE TUBE Ø6X10MM LENGTH 10M FOR ASPIRET, ASKIR



Ref. **BCASPI143**

### Description

Silicone tube roll diam 6 mm / ext 10mm x 10m

### Features

CE Marking	CE0123	Medical device	No
Bar code	0000000000888	HS Code	90330090
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	21	Packed unit / Height (cm)	5
Packed unit / Length (cm)	23	Packed unit / Weight	0.694
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	10	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

29/08/2024 17:35:59 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)