

FLOVAC SET 1L FOR ASKIR 118 ET 36BR



Ref. BCASPI501



29/08/2024 20:56:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

FLOVAC SET 1L FOR ASKIR 118 ET 36BR



Ref. BCASPI501

Description

Kit Flovac 1L, composed of: 1x BCASPI102: Jar Flovac 1L reusable for Aspiret, Askir 20 and 30 - 31843 1x BCASPI99: Sterile UU pocket 1L Flovac Aspiret, Askir 20 30 230 12V, 36BR 31848 1x

BCASPI135:
Set tubes Ø8x14mm +connectors Flovac ASKIR 36BR/C30/118/NEW HOSPIVAC 1x
BCASPI105: Round Retractor for System Flovac New Askir 20,30,230/12VBR
1x BCASPI164: Yellow Male Connector FLOVAC System (5 pcs) sf Askir12V/Emivac
Manual Configuration Aspiration Cleaners

Features

CE Marking	NONE	Medical device	No
Bar code	3661040191683	HS Code	4202121900
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.638
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 20:56:39 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr