

FLOVAC SET 2L FOR ASKIR 20/30/PROXIMITY/230/ASPIRET




Ref. BCASPI502



29/08/2024 20:57:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

FLOVAC SET 2L FOR ASKIR 20/30/PROXIMITY/230/ASPIRET



Ref. BCASPI502

Description

Kit Flovac 2 L, consisting of: 1x BCASPI103: Reusable Flovac 2L jar for Aspiret Askir 20, 30 and C30BR 31841 1x BCASPI100: Sterile UU pocket 2L Flovac Aspiret, Askir 20,30, 23012V,

C36BR 31846 1x

BCASPI111: Set tubes Ø6X10mm + connectors Flovac ASKIR 20/30/230/ASPIRET 1x

BCASPI105: Round Retractor for Flovac System New Askir 20,30,230/12VBR

1x Yellow Male Connector FLOVAC System (5 pcs) sf Askir12V/Emivac

1x Instructions for Configuration Suction Units

Features

CE Marking	NONE	Medical device	No
Bar code	3661040191690	HS Code	4202121900
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 20:57:32 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr