

FLOVAC SET 2L FOR ASKIR 118 ET 36BR



Ref. BCASPI503



29/08/2024 20:58:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

FLOVAC SET 2L FOR ASKIR 118 ET 36BR



Ref. BCASPI503

Description

Flovac 2L kit, consisting of: 1x BCASPI103: Reusable Flovac 2L jar for Aspiret Askir 20, 30 and C30BR 31841 1x BCASPI100: Sterile UU pocket 2L Flovac Aspiret, Askir 20,30, 23012V, C36BR 31846 1x

BCASPI135: Set tubes Ø8x14mm +connectors Flovac ASKIR 36BR/C30/118/NEW HOSPIVAC 1x

BCASPI105: Round Spreader for System Flovac New Askir 20,30,230/12VBR

1x Yellow Male Connector FLOVAC System (5 pcs) sf Askir12V/Emivac

1x Instructions for Configuration Suction Units

Features

CE Marking	NONE	Medical device	No
Bar code	3661040191706	HS Code	4202121900
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.844
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 20:58:37 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr