

ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L



Ref. BCFLOV10



29/08/2024 21:03:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L



Ref. BCFLOV10

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

Description

Ca-mi New Askir 20 Phlegm Vacuum - FLOVAC® Single Use Bottle/Bag, 1 L

Fixed electric surgical Suction Unit, designed to drain body fluids and perform oral, nasal and tracheal aspiration in adults and children. Suitable for tracheostomy patients and is ideal for minor surgical applications as well as postoperative home treatments.

Main applications: ENT, Home care, Room, Tracheostomy.

- Composed of: New Askir 20 Suction Unit + Flovac 1L Kit (1x Reusable Flovac 1L Jar / 1x Sterile UU pocket 1L Flovac / 1x Set tubes Ø6X10mm + connectors / 1x Round Retractor for Flovac System / 1x Yellow Male Connector FLOVAC System (5 pcs) / 1x Configuration Manual Suction Units

- The main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Optional roller foot

- Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug

- Engine: Oil-free piston pump requiring no maintenance

- Power supply: 230V-50Hz - other configurations available on request

- Maximum vacuum (adjustable): 0.75 Bar -75 kPa -563 mmHg

- Max free air flow: 16 l/min

- Noise level: 60.5 dBA

- Operating

cycle: Uninterrupted

- Weight: 2.5 kg / Dimensions: 35 x 21 x 18 cm

 29/08/2024 21:03:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L



Ref. BCFLOV10

Features

| | | | |
|------------------------------------|---------------|---------------------------|-----------|
| CE Marking | CE0123 | Medical device | Class IIa |
| Bar code | 3661040200019 | HS Code | 90189084 |
| Warranty (months) | 24 | Country | Italy |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | nc | Packed unit / Height (cm) | nc |
| Packed unit / Length (cm) | nc | Packed unit / Weight | 2.5 |
| Bare Unit / Weight (kg) | 2.5 | | |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 45 |
| Minimum humidity | 1 | Maximum humidity | 85 |
| OTHER FEATURES | | | |

29/08/2024 21:03:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

ASKIR CA-MINEW SUCTION UNIT 20 BOTTLE/BAG UU FLOVAC 1L



Ref. BCFLOV10

Content

| Ref | Designation | Quantity |
|-----------|---|----------|
| BCASPI99 | DISPOSABLE FLOVAC® LINER, 1000 ML for Aspiret, Askir 20, 30 | 1 |
| BCASPI164 | Male Connector (for Vacuum Port) for Askir 20 pack of 5's | 1 |
| BCASPI102 | Container Flovac 1L (reusable) for Aspiret, Askir 20 and 30 | 1 |
| BCASPI504 | IFU for CA-MI suction units | 1 |
| BCASPI06 | New Askir 20 electrical suction unit, 1L autoclavable polycarbonate j | 1 |
| BCASPI111 | Set of Silicone Tubes Ø 6x10 Flovac + connector Aspiret, Askir 20/30 | 1 |

29/08/2024 21:03:36 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr