

## CA-MI NEW ASKIR SUCTION UNIT 30 BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV17



29/08/2024 21:05:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW ASKIR SUCTION UNIT 30 BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV17

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

### Description

Capom New Askir 30 - FLOVAC® Single Use Bottle/Bag, 2000 ml

Fixed electric surgical Suction Unit, intended to drain body fluids and to perform oral, nasal and tracheal aspiration in adults and children.

Suitable for tracheostomy patients and is ideal for minor surgical applications as well as postoperative home treatments.

Main applications: Endoscopy, Gastroenterology, General surgery, Home care, Obstetrics, Room, Veterinary Med.

Composed of: New Askir 30 Suction Unit + Flovac 2L Kit (1x Reusable Flovac 2L Jar / 1x Sterile UU pocket 2L Flovac / 1x Set tubes Ø6X10mm + connectors / 1x Round Retractor for Flovac System / 1x Yellow Male Connector FLOVAC System (5 pcs) / 1x Configuration Manual Suction Units

The main system is equipped with a pressure gauge (bar & kPa) and a vacuum regulator

- Optional roller foot

- Supplied accessories: Fluid collection jar with an anti-reflux valve with different options / 6 x 10 Silicone hoses (autoclavable) and conical connectors ø 8-9-10 mm / Antibacterial & hydrophobic filter (single patient) / CH20 cannula / Power cord with Schuko plug

Motor: Oil-free piston pump requiring no maintenance

Power supply: 230V-50Hz - other configurations available on request

- Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

Max free air flow: 40 l/min

Noise level: 60.5 dBA

- Uninterrupted operating cycle

Weight: 3.6 kg / Dimensions: 35 x 21 x 18 cm

29/08/2024 21:55:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## CA-MI NEW ASKIR SUCTION UNIT 30 BOTTLE/BAG UU FLOVAC 2L



Ref. BCFLOV17

### Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	3661040200026	HS Code	90189084
Warranty (months)	24	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	3.6
Bare Unit / Weight (kg)	3.6		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			

29/08/2024 21:05:35 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)