


## DERMAL CURETTES 4MM BOX OF 10

Ref. CPCURE03



29/08/2024 20:30:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## DERMAL CURETTES 4MM BOX OF 10



Ref. CPCURE03

Medical device used in light skin surgery procedures

### Description

The very tapered sharp edge ring allows tissue harvesting and provides good visibility of the operating field.

The ergonomic handle allows a good grip of the instrument and ensures a safe and precise gesture.

- Equipped with a

- cap to prevent any risk of cut

- Sterilization by gamma irradiation 2.5 MEGA-RADS

- Packaging in individual sterile peelable bag covered with a transparent film

- Blue sleeve

- Length of the handle: 13 cm

- Box of 10

- Ready to use

- Complies with European directives 93/42/EC

### Features

CE Marking	CE2460	Medical device	Class IIa
Bar code	3661040186016	HS Code	90183290
Warranty (months)	12	Sterile	Yes
Country	India		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	4	Packed unit / Height (cm)	7
Packed unit / Length (cm)	21	Bare Unit / Length (cm)	13.4
Packed unit / Weight	0.144		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	10	Max temperature	40
Minimum humidity	35	Maximum humidity	65

#### OTHER FEATURES

29/08/2024 20:30:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)