


## PUNCH BIOPSY Ø 3,5 MM BOX OF 10

Ref. CPPUNC06



29/08/2024 20:15:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

[www.gsh-med.fr](http://www.gsh-med.fr)

+33(0)4 42 90 31 31

[info@gsh-med.fr](mailto:info@gsh-med.fr)

## PUNCH BIOPSY Ø 3,5 MM BOX OF 10



Ref. CPPUNC06

Ready-to-use, single-use device

### Description

Efficient and precise thanks to its very sharp seamless stainless steel border. The notched neck ensures a good grip of the instrument. The size of the punch indicated on the handle allows easy identification. - Dimensions of the

handle: L 82 x diameter 10 mm

- Sterilization by gamma irradiation higher than the standard of 2.5 MEGA-RADS

- Under individual sterile peelable packaging covered with a transparent film

- Box of 10

- Complies with European directives 93/42/EC

### Features

CE Marking	CE2460	Medical device	Class IIa
Bar code	3661040186207	HS Code	90183290
Warranty (months)	12	Sterile	Yes
Country	India		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Bare Unit / Length (cm)	8.2
Packed unit / Weight	0.064		
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	10	Max temperature	40
Minimum humidity	35	Maximum humidity	65

### OTHER FEATURES

📅 29/08/2024 20:15:20 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)