

AUTOCLAVABLE NYLON BRUSH (12 PIECES)

Ref. DB01201_CCC



 29/08/2024 21:43:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

AUTOCLAVABLE NYLON BRUSH (12 PIECES)



Ref. **DB01201_CCC**

Description

- Surgical brush
- Gentle but effective cleaning of surgical instruments and nails
- Nylon brush bristles, very resistant to remove visible dirt (traces of dried blood, fabrics, rust film ...).
- Autoclave sterilization at 100°C.

Features

CE Marking	NONE	Medical device	No
Bar code	8023279257151	HS Code	96032980
Warranty (months)	12	Country	China
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.55
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 21:43:04 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr