

THERMAL PRINTER FOR KERN- BALANCES WITH DATA INTERFACE RS-232

KERN®


Ref. EKMWA07



 29/08/2024 17:30:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

THERMAL PRINTER FOR KERN-BALANCES WITH DATA INTERFACE RS-232



Ref. **EKMWA07**

Description

- Date and time printing
- Thermal printing

- Paper strips: Width 112 mm
 - Case dimensions L×P×H 165×140×50 mm
 Tip: buy the scale and printer together

Features

CE Marking	NONE	Medical device	No
Bar code	3661040173702	HS Code	84239010
Warranty (months)	12	Country	Germany
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	8
Packed unit / Length (cm)	30	Packed unit / Weight	1
Bare Unit / Weight (kg)	0.35		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	NA	Max temperature	NA
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:30:49 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
 13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr