

## PLASTIC TROLLEY, 2 TRAYS 61,5 X 45 CM


Ref. OM614501



 29/08/2024 21:14:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## PLASTIC TROLLEY, 2 TRAYS 61,5 X 45 CM

Ref. **OM614501**

**Lightweight, handy and very robust pedestal table**

### Description

The quality of the casters and ergonomic handles ensure maximum comfort when travelling.

External dimensions: 80 x 48 x H 95 cm, including handles

- Dimensions of the trays: 61.5 x 45 cm

-Maximum load: 180 kg per trolley and 60 kg per level

-Color: Light gray

-Materials: Aluminum upright, PP top

- 4 PVC swivel casters, Ø 100 mm including 2 brake

-Weight: 8.5 kg


### Features

CE Marking	NONE	Medical device	No
Bar code	3760206421227	HS Code	8716800000
Warranty (months)	12		
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	71	Bare Unit / Width (cm)	48
Packed unit / Height (cm)	19	Bare Unit / Height (cm)	95
Packed unit / Length (cm)	52	Bare Unit / Length (cm)	80
Packed unit / Weight	9.806	Bare Unit / Weight (kg)	8.5
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	5	Max temperature	35
Minimum humidity	NA	Maximum humidity	NA
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

 29/08/2024 21:14:30 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)