

PLASTIC TROLLEY, 3 TRAYS 61,5 X 45 CM - WITH HANDLES

Ref. OM614502




CE



 29/08/2024 21:14:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

PLASTIC TROLLEY, 3 TRAYS 61,5 X 45 CM - WITH HANDLES

Ref. **OM614502**

Lightweight, handy and very robust pedestal table

Description

The quality of the casters and ergonomic handles ensure maximum comfort when travelling. - External dimensions: 80 x 48 x H 100 cm,

including handles

- Dimensions of the trays: 61.5 x 45 cm
- Maximum load: 180 kg per trolley and 60 kg per level
- Color: Light gray
- Materials: Aluminum upright, PP top
- 4 PVC swivel casters, Ø 100 mm including 2 brake
- Weight: 8.5 kg

Features

CE Marking	CE	Medical device	No
Bar code	3760206421234	HS Code	8716800000
Warranty (months)	12		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	52	Bare Unit / Width (cm)	48
Packed unit / Height (cm)	30	Bare Unit / Height (cm)	100
Packed unit / Length (cm)	78	Bare Unit / Length (cm)	80
Packed unit / Weight	13	Bare Unit / Weight (kg)	8.5
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	5	Max temperature	35
Minimum humidity	NA	Maximum humidity	NA
OTHER FEATURES			
Packed unit / Quantity	1		

 29/08/2024 21:14:57 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr