

**5000ML AUTOCLAVABLE
POLYCARBONATE JAR ASKIR NEW
ASPIRET NEW HOSPIVAC**




Ref. BCASPI108

 29/08/2024 17:33:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

5000ML AUTOCLAVABLE POLYCARBONATE JAR ASKIR NEW ASPIRET NEW HOSPIVAC



Ref. **BCASPI108**

Features

CE Marking	NONE	Medical device	No
Bar code	0000000000659	HS Code	4202121900
Country	Italy		
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	nc
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:33:53 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr