

## SET OF SILICONE TUBES Ø8X14 JAR + CONNECTOR FOR ASKIR C30BR



Ref. BCASPI115



29/08/2024 15:52:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

 [www.gsh-med.fr](http://www.gsh-med.fr)

 +33(0)4 42 90 31 31

 [info@gsh-med.fr](mailto:info@gsh-med.fr)

## SET OF SILICONE TUBES Ø8X14 JAR + CONNECTOR FOR ASKIR C30BR



Ref. BCASPI115

### Description

Set of pipes containing:- 1 pipe x diam 8 / ext 14mm x length

150

cm-1 pipe x diam 8 / ext 14mm x length 34 cm

Conical connectors 10-11-12 mm

### Features

CE Marking	CE0123	Medical device	No
Bar code	3661040151069	HS Code	90330090
Warranty (months)	12	Country	Italy
<b>DIMENSIONS AND WEIGHT</b>			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	0.256
<b>SPECIFIC STORAGE CONDITIONS</b>			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
<b>OTHER FEATURES</b>			
Packed unit / Quantity	1		

📅 29/08/2024 15:52:55 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand  
13290 Aix-en-Provence | France

🌐 [www.gsh-med.fr](http://www.gsh-med.fr)

☎ +33(0)4 42 90 31 31

✉ [info@gsh-med.fr](mailto:info@gsh-med.fr)