

LID FOR 5000ML POLYCARBONATE JAR NEW HOSPIVAC



Ref. BCASPI120



 29/08/2024 17:34:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

LID FOR 5000ML POLYCARBONATE JAR NEW HOSPIVAC



Ref. BCASPI120

Features

CE Marking	NONE	Medical device	No
Bar code	0000000000710	HS Code	4202121900
Warranty (months)	12	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	18	Packed unit / Height (cm)	12
Packed unit / Length (cm)	18	Packed unit / Weight	0.272
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			
Packed unit / Quantity	1		

29/08/2024 17:34:18 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr