

**ROLL OF SILICONE TUBE
Ø8X14MMX50M ASKIR
36BR/118/C30/C30BR/N-HOSPIVAC**



Ref. BCASPI146



 29/08/2024 17:38:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

ROLL OF SILICONE TUBE Ø8X14MMX50M ASKIR 36BR/118/C30/C30BR/N-HOSPIVAC



Ref. BCASPI146

Description

Roll silicone tube diam 8mm / ext 14mm x 50 m

Features

| | | | |
|------------------------------------|---------------|---------------------------|----------|
| CE Marking | CE0123 | Medical device | No |
| Bar code | 0000000000901 | HS Code | 90330090 |
| Warranty (months) | 12 | Country | Italy |
| DIMENSIONS AND WEIGHT | | | |
| Packed unit / Width (cm) | 42 | Packed unit / Height (cm) | 14 |
| Packed unit / Length (cm) | 42 | Packed unit / Weight | 3.45 |
| SPECIFIC STORAGE CONDITIONS | | | |
| Minimum temperature | -20 | Max temperature | 45 |
| Minimum humidity | NA | Maximum humidity | 85 |
| OTHER FEATURES | | | |
| Packed unit / Quantity | 1 | | |

29/08/2024 17:38:24 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr