

SUCTION UNIT CA-MI ASKIR 36 BR POLYCARBONATE JAR 1L



Ref. BCASPI24



29/08/2024 14:47:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

www.gsh-med.fr

+33(0)4 42 90 31 31

info@gsh-med.fr

SUCTION UNIT CA-MI ASKIR 36 BR POLYCARBONATE JAR 1L



Ref. **BCASPI24**

Designed for professional use and equipped with a powerful suction system that can drain up to 36 l/min max, Askir 36BR is lightweight and easy to carry.

Description

Suction Unit Ca-mi Askir 36 BR, polycarbonate jar 1000 ml

, 36 l / min

Designed for professional use and equipped with a powerful suction system that can drain up to 36 l/min max, Askir 36BR is lightweight and easy to carry.

Main applications: Ambulance car, Emergency Department, Home care, patient transfer, Room, Tracheostomy, Veterinary Med.

- Three different operating options: The AC/DC adapter, rechargeable battery and 12V unit make this surgical Suction Unit a versatile unit guaranteeing CONTINUOUS interventions without the risk of overheating.

- Equipped with a visual and audible alarm system indicating if the battery level is low, the main system is also equipped with a pressure gauge (bar & kPa) and a vacuum regulator.

-

Foot with wheels and optional carrying case

-

Accessories provided: Makrolon 2858 autoclave polycarbonate fluid collection jar, equipped with an anti-reflux valve with different options / 8 x 14 Silicone hoses (autoclavable) and conical connectors ø 10-11-12 mm / Antibacterial & Hydrophobic filter (single patient) / Cannula CH20-

Power cord with Schuko plug, Automotive Adapter 12V,

Universal adapter AC / DC-

Motor: Oil-free piston pump requiring no maintenance

- Power supply: 14V 4A with Universal AC / DC adapter / rechargeable Internal battery (12V 4A) Autonomy 60min

- Charging time 240min / Automotive adapter (12V 4A

)-

Maximum vacuum (adjustable): 0.80 Bar -80 kPa -600 mmHg

-

Max free air flow: 36 l/min (without jar)

- Noise level: 63 dBA

 29/08/2024 14:47:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr

SUCTION UNIT CA-MI ASKIR 36 BR POLYCARBONATE JAR 1L



Ref. BCASPI24


Features

CE Marking	CE0123	Medical device	Class IIa
Bar code	8054619100424	HS Code	90189084
Warranty (months)	36	Country	Italy
DIMENSIONS AND WEIGHT			
Packed unit / Width (cm)	nc	Packed unit / Height (cm)	nc
Packed unit / Length (cm)	nc	Packed unit / Weight	5.5
Bare Unit / Weight (kg)	4.39		
SPECIFIC STORAGE CONDITIONS			
Minimum temperature	-20	Max temperature	45
Minimum humidity	NA	Maximum humidity	85
OTHER FEATURES			

 29/08/2024 14:47:21 - For information purposes only, can be subject to change

GSH | 30 rue Jean de Guiramand
13290 Aix-en-Provence | France

 www.gsh-med.fr

 +33(0)4 42 90 31 31

 info@gsh-med.fr